



P.O. Box 264, Westfield, Indiana 46074

Contact your local registered salesperson who contacted you first for more information, or

Call: 317-804-1803 Email: info@tlchealthsystems.com Website: <http://www.tlchealthsystems.com>

Pain Relief and other Health Benefits at the Speed of Light – Patent Pending – Scientifically Engineered

Owner's Manual for TLC Spectrabrite™ Focused Light Disks and Pens

Treatment Comparison

If there were a way you could own the TLC Spectrabrite™ Disk or Pen and use it as often as desired for pain relief (without side effects!) for the cash price of one month's supply of prescription painkillers or pain patches (with their stated potentially-*dangerous* side effects), which one would you choose?



Buy a One-Month's Supply of Painkillers and/or Pain Patches



Or buy a TLC Spectrabrite™ Disk – to use for years



And also own a TLC Spectrabrite™ Pen to use for years to come



P.O. Box 264, Westfield, Indiana 46074

Contact your local registered salesperson who contacted you first for more information, or
Call: 317-804-1803 Email: info@tlchealthsystems.com Website: <http://www.tlchealthsystems.com>

Pain Relief and other Health Benefits at the Speed of Light – Patent Pending – Scientifically Engineered

Dear Valued Customer,

It is an honor to provide and serve you with our amazing and very unique TLC Spectrabrite™ technology that has taken over 20 years of research to perfect. We want to help you realize its potential for you, your family, friends, clients, patients, and even your pets!

Each of our TLC Spectrabrite™ products has been designed and engineered with many proprietary health programs. Each program is delivered via: light pulsing frequencies, sine waves, longitudinal scalar wave energy, and sound, whether you use the preprogrammed settings or manually set colors to address certain issues. (Instructions provided in the manual)

Light is an essential food for plants; likewise, we too, require light for health. Now by using TLC Spectrabrite's™ technology you may access all the beneficial light frequencies while eliminating harmful rays and bad frequencies. This allows you to enjoy almost innumerable TLC Spectrabrite's™ light therapies that are beneficial to your health, your family, friends and patients. So many of our customers have claimed TLC Spectrabrite™ has helped relieve their patients from sicknesses and pains without side effects. We at TLC desire this for you and everyone else as well!

TLC Health Systems uses Spectrabrite only as our delivery system for certain amounts of our TLC Health Systems scientific and proprietary health protocols. If what you are about to purchase is just a Spectrabrite look alike from someone else and it is not our official TLC Health System "TLC Spectrabrite™" product(s), then you are not getting the real thing. Any substitute will not do what our complete technology does. We alone have what is called TLC Spectrabrite™ products and no other look alike imitations even with the Spectrabrite's name on it can compare with or do what our official TLC Spectrabrite™ products do.

TLC Health Systems also specializes in Class III and Class IV lasers as well.

TLC Spectrabrite™: Pain Relief and other Health Benefits at the Speed of Light.

From all of us at,

T.L.C. Health Systems
TLC Health Systems



Pain Relief and other Health Benefits at the Speed of Light – Patent Pending – Scientifically Engineered

Owner's Manual - For Your Unique TLC Spectrabrite™ Pens and Disks

Letter From Our TLC Spectrabrite™ Inventor

The use of light therapy is a very old concept indeed. Many ancient cultures practiced various forms of Helio-therapy, including the people of ancient Egypt, Greece and Rome. The Inca, Assyrian and early Germanic tribes also used the sun for health-augmenting purposes. Also, Indian medical literature dating back to 1500 BC describes several treatments which combine herbs and natural sunlight to treat assorted skin conditions. Similar references are also found in Oriental writings dating back to 200AD.

However, a Faroese physician, one Niels Finsen, is believed to be the father of modern phototherapy. He developed the first artificial light source for this purpose, and used his invention to treat lupus vulgaris. He was awarded the Nobel Prize in Physiology in the year 1903. Since that time, large arrays of treatments have been developed from the use of controlled light sources, and various light frequencies. These devices are commonly used to treat such ailments as acne and swelling, by means of low-level laser light, red light, near infrared and ultraviolet frequencies. However, relief from pain constitutes the greatest health utilization of this technology.

Most apparatuses designed for light therapy applications consist of a single light source, emitting a single color, and a single frequency band. Most of these light sources emit steady-state radiation and are designed for low energy application.

The present version of our TLC Spectrabrite™ Clinical Unit totally eclipses all existing technologies in the lower price ranges, and supplies greater energy delivery, selectable programmed frequencies, and the versatility of multiple wavelength utilization, by program choice, wand selection, or both. The apparatus also contains a micro-processor which is capable of holding 999 light therapy programs. These programs automatically modify the behavior of the entire instrument so that it is capable of performing a multitude of tasks, such as changing basic color tones, altering pulse rates, blending colors automatically, timing each light therapy session, and much more.

We at Millennial Health Systems LLC along with our partners at TLC Health Systems have upgraded and improved the TLC Spectrabrite™ units to be the most advanced units of their kind and because there are no other Spectrabrite devices like the unique TLC Spectrabrite™ units. We are confident you will benefit greatly from using any of their Unique TLC Spectrabrite™ products, whether it is the TLC Spectrabrite™ Clinical Unit, their Light Pens, Light Disks, etc. All can deliver years of trouble free service, and are reasonably simple to use. We wish you good health and blessings on your journey "Into the Light."

Welcome to the TLC Health Systems Family of TLC Spectrabrite™ units containing Focused Low Level Light Therapy (LLLT) Instruments.

Inventor,

David R. High

P.S. Please read this entire manual before using your unique TLC Spectrabrite™ equipment.



P.O. Box 264, Westfield, Indiana 46074

Contact your local registered salesperson who contacted you first for more information, or
Call: 317-804-1803 Email: info@tlchealthsystems.com Website: <http://www.tlchealthsystems.com>

Pain Relief and other Health Benefits at the Speed of Light – Patent Pending – Scientifically Engineered

Owner's Manual for TLC Spectrabrite™ Disk and Pen

BASIC HISTORY: Our TLC Spectrabrite™ Focused - Low Level Light Therapy – **Disk** was made originally as a toy and our Focused - Low Level Light Therapy – **Pen** was initially made for pain care. But our associates, doctors, researchers, inventors and scientists wanted to make them both to carry our TLC Spectrabrite™ special proprietary health technology, so we re-engineered everything. We scientifically programed all our proprietary info, patent pending processes, and exclusive signals that are recognized for reducing inflammation, pain and much more into our new Focused – Low Level Light Therapy (LLLT) devices and they work very well.

TECHNOLOGY - HOW IT WORKS: This programming process involved infusing our laser pump diode, vanadate crystal, dielectric mirror, copper heat sink, potassium crystal, cavity resonator, and collimator lens with longitudinal scalar wave energy and information. The longitudinal scalar wave segment of the wave equation is not in current textbooks but was developed and tested over the years by NASA, many Medical Schools, scientists, doctors and researchers worldwide, including but not limited to, Dr. William Campbell Douglass, MD author of “Into the Light”, Darius Dinshah, S-CN author of “Let There Be Light”, Nikola Tesla, and many more found on unbiased websites like American Institute of Light Therapy www.aiolt.org along with other researchers and doctors who will explain more than we can here, but what we do write tells our story.

LIGHT BEAMS and COLORS: All TLC Health Systems devices utilize beams of light and color in our Focused Light Therapy Disks and Pens, commonly referred to as *Low-level Laser Therapy (LLLT)*, *Laser Therapy*, *LED Focused Light Therapy*, and *Phototherapy* to deliver our therapeutic protocols and knowledge to injured areas. These LLLT technologies assist us in delivering our exclusive, patent pending, proprietary protocols, that no one else in the world has, and they do so by transporting our longitudinal scalar wave protocols, like Nikola Tesla and many other well-known scientists discovered, to help deliver our Dunamis Power to damaged areas, offering pain relief at the speed of light, and with diligent use of our LLLT you may find other beneficial situations occurring, helping even more.

LONGITUDINAL SCALAR WAVES: The key to longitudinal scalar waves is they propagate similar to a plasma wave, acoustic wave or gravitational wave and have the ability to penetrate any material including a Faraday Cage. The longitudinal wave consists of particles or formations that can carry information (impulses) and energy in the same wave form. This process is carried out in a transponder system in which a signal can be directed in a controlled way into material that is subject to programming holograms, silicone, crystals, liquids and other materials that can hold a memory through our Focused Light Therapy Device(s) and many of our other units too, especially “TLC Spectrabrite™”.

LIGHT CONVERTED TO SOUND and SOUND BACK TO LIGHT: By understanding that light and sound correlate and that light can be converted to sound and sound to light gives the opportunity to take (nm) of light, convert to sound in (Hz), further converting to longitudinal scalar waves and imprint them into the material that makes up our LLLT. Every molecule has an emission signature of light in nano meters (nm), so we can reproduce the molecule as light taking it directly to the body's area of need, so it will react to the light beam in this case through the combination of signals being beamed through our TLC Spectrabrite™ units. We then just have to select the molecules that deal with pain and inflammation and imprint into the material that makes up our LLLT and the body responds to these focused light beams that carry our protocols. Our LLLT use multiple light beams to carry and deliver Pain Therapy Tesla, Pain Therapy 528, Spectra Schumann and dozens of other proprietary therapies that electronically pulse deep and radiate out to reach painful and/or injured areas. Light Pen and Disk's Output is 5 mW and wavelength 650 nm or less.

HOW TO USE: Your LLLT unit uses multiple colors of light to carry our proprietary health frequencies.



TLC Health Systems

1. Red: Pain Relief
2. Green: Healing Properties
3. Blue: Youth and Beauty - Acne – Repair – Calming – Blood-borne diseases – Rheumatoid Arthritis
4. Gold: Helps dissolve blood clots, and produce digestive fluids
5. Yellow: Energy, Lymphatic system
6. Lime: Bone Builder, Nutrition
7. Aqua: "Skinny" Appetite Control
8. Sky: Brain Boost, Memory and Focus
9. Cyan: Stress Relief, Peace and Calm
10. Purple: Sleep, Good Dreams
11. Violet: Passion – Spleen Builder
12. Magenta: Happy Blissful Sleep Improvement – opens blood vessels
13. White: Healing - wellbeing overall
14. Fade: slowly color to color on auto
15. Jump: from color to color on auto
16. Reset: Start from the beginning
17. Decreases Light Intensity
18. Increases Light Intensity

TLC Spectrabrite™ Focused Light Therapy Devices (LLLT) are used by Doctors, including M.D., Chiropractors, Dentists, Veterinarians and others. By using our low-level light-emitting diodes to help alter cellular function it brings welcome pain relief, wellness and there are no known side effects. **Disk Batteries** - 3 (AAA) and Remote - 2 (AS 110 coin batteries). Disk AAA batteries may last up to 13 hrs. under constant use but even if lights are still lit remote may not turn off Disk so change 3 AAA Disk batteries only and NOT the remote batteries. If used often, buy recharge-able batteries and a charger. **CAUTION:** To prevent overheating turn pens or disks off at least every 2 hours, let them cool down, and then restart.

FROM THE LAYMAN'S PERSPECTIVE: Use your TLC Spectrabrite™ Focused Light Therapy Device (LLLT) **on any area of pain as often and as long as you need to.** You can start by using the LLLT in a circular or sweeping motion above the skin over the general area of pain, or directly on the skin for several minutes. This may help to reduce inflammation. Then massage the area with your fingertips to locate the most painful points. Put the Disk directly on that skin area after identifying those points (you can also use your LLLT directly on the inflamed or fresh wound areas). Hold LLLT on the pain area, meridian, acupuncture or reflexology point, or any other point for five minutes or five hours if you wish. **There is no maximum amount of time that you may spend on any painful point.** If you feel any responses from the body such as tingling, warmth, muscle responses, a sensation going down a limb or a response in another area of your body, work on that point until the responses subside. This is an indication that you are on a point that needs the attention of your LLLT. **You cannot use your LLLT too much or too often.** Study impartial websites www.aiolt.org save to favorites. Search LLLT and sickness name(s) as to whether or not LLLT technology worked effectively or not on whatever the sickness or symptoms name may be that needs tending. **CAUTION:** To avoid overheating, turn pens or disks off at least every 1 - 2 hours if necessary, let units cool down, and then restart.

SPINE APPLICATION: Apply LLLT directly on the skin approximately 1 ½ inches from the spine on either side of the first two vertebrae (C1 and C2). Doing this may produce more beneficial results and also in the area of the “meridians” which are of greater electrical conductivity in the body and thus more good things may occur.

One of the most important areas of the body to use your LLLT on other than the areas of primary pain is at the base of the skull in the area of the brainstem and all the way down the spine to the tailbone and back up both sides of the spine again. The first two vertebrae (C1 and C2) influence the entire body and the brainstem is the control center for the central nervous system. Reducing inflammation in this area may be quite beneficial. Hold Light Disk in place of pain, on spinal column, and for overall wellbeing directly on wrist where pulse is taken and do so for hours by wrapping it in place with an ACE bandage or Velcro strap. For more info on spine areas where LLLT has been applied positively influencing internal organs, view: American Institute of Light Therapy www.aiolt.org or: U.S. National Library of Medicine, or National Institute of Health: <http://www.nih.gov> or: https://en.wikipedia.org/wiki/Low_level_laser_therapy or: **(do online searches for LLLT and any sickness name)**. More uses for LLLT and other TLC patent pending items are discovered daily.

- ***Joint problems** such as injuries, arthritis, bad knees and other joints apply LLLT around the entire joint area.
- *Use your LLLT on **old injuries** such as **broken bones, sprains and scars.**
- *Use your LLLT on **warts, moles and wrinkles** (ten to twenty minutes every day).
- *Use LLLT to **whiten teeth** by shining directly on the teeth for ten minutes or more one or more times daily.
- *Use your LLLT on any skin condition such as **acne, psoriasis, shingles, rashes,** etc.
- *Use LLLT on **(tinnitus)** ringing in the ear, place LLLT directly over your ears for at least ten minutes a day.
- *Use LLLT on **nail fungus** for at least ten minutes one or more times per day, or fasten on with an ace bandage.
- *Use LLLT on **Headaches and Migraines** hold LLLT on the neck area between it and the ear for 30 minutes.
- *Use LLLT on all body parts issues and it is **recommended to be used on all food and beverages too.**

Your LLLT is antifungal, antibacterial and antiviral. For sure, use your LLLT on fruit juice, coffee, tea and other beverages (except carbonated drinks because it will cause them to go flat) for five minutes for a much smoother taste and to add many therapeutic values to every drink, just do it. For 5+ minutes set your (see thru) clear glass or cup with drink in it on top of your LLLT Disk; turn to jump, fade or any color desired, then drink.

MORNING - DRINKING WATER EXERCISE:

Shine your LLLT on your drinking water every morning for 5 minutes. Before eating or drinking anything in the morning shine LLLT on 24 ounces of water for five minute or more on any desired setting. If you have powdered or liquid minerals and vitamins add them to your water, then LLLT the water from the top, or if you have a (see through glass vessel) set your drink on your LLLT **Disk** and turn unit to any setting. By adding powdered or liquid minerals and vitamins to your water you'll find it will help to retain the health signals programmed into your LLLT longer which are carried into our body's blood cells along with much needed

water, vitamins and minerals which helps to prevent being dehydrated. Hydration of our cells in the morning is very important. If possible don't eat or drink anything for 30 to 45 minutes after drinking your LLLT water.

By doing this morning water exercise faithfully and throughout day, it may help marvelously reduce internal and external pains and increase your energy level throughout the day and much more. If you can't drink that much water, begin with as much as you can drink and increase as you go. To get the maximum benefit from your LLLT the body must be hydrated and our mineral and vitamin supplies adequate, thus by LLLT-ing our water each morning and even throughout the day may help you immensely! **We recommend our TLC Ionic, Fulvic, Humic Minerals, Vitamins and Amino Acids** 1 oz. per glass of LLLT water, there's nothing else like it!

TREATMENTS: Doctors and Chiropractors after using our Spectrabrite Clinical Unit on their patients while in their offices, will then treat with the Focused Low Level Light Therapy Disk and/or both Focused Light Therapy Pen for their patients to rent and use faithfully in the privacy of their own homes until their next office visit.

A Few Conditions Where Low-Level Laser (or Light) Therapy (LLLT) Are Therapeutically Applied By Doctors

GENERAL

Repetitive Stress Injury
Carpal Tunnel Syndrome
Epicondylitis, Fibromyalgia
Reflex Sympathetic Dystrophy (RSD)
Rotator Cuff Tear, ACL Tears
Temporo-Mandibular Joint Dysfunction (TMJ)

INJURIES

Ligament and strained Tendon
Fractures with Associated Soft Tissue Injuries
Facet Joint Syndrome
Contusions
Bulging and Herniated Discs – Broken Bones

INFLAMMATORY

Tendonitis, Bursitis, Synovitis
Bad Knees, Hips, Shoulders, and other Joints
Myositis - Autoimmune conditions
Plantar Fascitis
Rheumatoid Arthritis, any type Arthritis / Gout
Shingles (herpes zoster)
Eczemas, Psoriasis
Teeth - pain, bleeding, gums gingivitis and etc.

Acne, Acne Rosacea
Genital Herpes, viruses, pathogens, parasites, fungus
MS – Multiple Sclerosis

OTHER APPLICATIONS

Eyes - Floaters may vanish w/applied on feet and wrists
Back – Vertebrae's (C1 and C2) and entire spinal column
Wound Healing (acute and chronic)
Scars – even from decades before
Dermal Ulcers
Venous Stasis
Atherosclerotic
Contact Diseases
Diabetic
Traumatic Injuries
Burns: Thermal, Chemical
Headaches, Migraines
Lymphedema
Neuropathies (diabetic, neurogenic)
Addictions
Warts, skin rashes
Remodeling Dermal Lesions (keloid scars) and more!

Search (LLLT and sickness names) find unbiased sites: www.AIOLT.org

Many researchers use LLLT on Foot and Hand Reflexology points and on Acupuncture points as well.

WARRANTY: All TLC LLLT carry a Ninety (90) Day Warranty from the date of delivery against any defects. ***Do not drop your TLC LLLT on the ground, floor, or in water even though it's designed to be waterproof.***

BATTERIES: Unscrew Light Disk to replace 3 (AAA) and your Remote requires two (AS 110 coin batteries).

TROUBLE SHOOTING: If Disk lights are still working, meaning lit up, but your remote won't turn the lights off, then it is probably NOT the remote batteries **BUT** the three AAA batteries that need replacing 99% of the time. So, first try replacing the three AAA batteries in Disk and not the two Remote (AS 110 coin batteries). Rarely do the remote control batteries need replacing, but if they do, replace with two (AS 110 coin batteries).

NOTICE: As wonderful as our TLC Spectrabrite™ LLLT devices are none of them can replace or compare with our more sophisticated, TLC Clinical Model used by doctors', hospitals, clinics, chiropractors, dentists, veterinarians and/or their patients in the privacy of their own home. To learn more about our TLC Spectrabrite™ units contact your representative or call us and we'll put you in touch with one in your area.

SPECTRABRITE'S CLINICAL UNIT IS THE ULTIMATE

Doctors, own your own ultimate TLC Health Care Systems unit. *Nothing else we know of can replace our TLC Spectrabrite™ Clinical Model with 4 rectangular wellness pads shown below, all 8 wands if you desire, for each wand has a specific purpose, but for sure get these 3 wands **RGB**, **UV** and **Amber** along with the 4 wellness pads shown below so you, your family, friends, and patients can be on the way to a more healthy life.*




Learn by studying Spectrabrite www.tlchealthsystems.com take the video tour, 3 min video, and/or do 3 hour doctors training.

Contact Your Local TLC Representative to Learn More about TLC Products or Contact Company

TLC Health Systems
P.O. Box 264 Westfield, IN 46074
Phone: 317-804-1803 Email: info@tlchealthsystems.com
Website: www.tlchealthsystems.com

TLC Health Systems
P.O. Box 264
Westfield, IN. 46074
317-804-1803
info@tlchealthsystems.com
www.tlchealthsystems.com



This picture is the Ultimate Package for Every Doctors' Office. All these other products shown are treatments for Their Valued Patients to utilize at home until their next scheduled Doctors' visit. There are no harmful side effects!

Just A Few Testimonies from People Who Have Used Our TLC Spectrabrite™ by TLC Health Systems

These testimonials below and many more will be signed and witnessed and held on file at our TLC home office

Diabetic – My feet swell up and hurt very badly and I cannot walk if I eat too much salt. One night we're out for dinner the food contained too much salt and by the time I got home I could hardly walk and in pain. I put the Disk on one foot before bed time and on the other foot while I slept. I woke up totally healed the next morning, no pain or swelling. *Lance*

Hernia – I was in much pain, I used the Disk everyday possible on my belly, holding disk in place with an Ace Bandage and my batteries lasted about 12 hours before replacing. Everything is back in place somehow and No more pain. *Karen T.*

Pneumonia – I get very sick at least once each year, coughing, fever, sore throat, aches and pains to the point of being bed ridden with Pneumonia, I know these symptoms all too well. I was getting worse, fever, lungs filling us, etc. so I put the disk on my throat fastened with an Ace Bandage and went to sleep. I woke up totally healed the next morning. My whole family at one time or another is using our Disk, and now we have a Pen too, almost 24 hours a day now. *LJ*

Back Pain – For 3 ½ years I have suffered with back pain with no relief. I never took pain killers because I did not want to get hooked, so I toughed it out to no avail, until one day a man said here put this disk on your back where the pain is, lean back in your chair while you work, 20 minutes later my pain was gone. I strap the disk on my body in different places and sleep with it on all night, and honestly I have not felt this good for years, it is wonderful. *Carolyn L.*

M.S. Multiple Sclerosis – My friend's sister, she suffered in pain with M.S. for years, unable to standup and walk. So, I told her to use this TLC-Pen on her sister. She did on her back for 1 hour the first day, starting from the top of her spin down to the tail bone and back up again and then back down both sides of the spine and she did so each day thereafter for an hour or longer if possible. Her sister literally stands up on her own now and walks with assistance. *Dave W.*

Arthritis – My hands have become so stiff over time that I can hardly use them and I have constant pain. A friend said hey, hold this disk in your hand and see if your pain goes away, and I said for how long, he said I don't know, I guess until the pain goes away. 20 – 30 minutes later my pain was gone so I switched hands. I use it every day for preventative maintenance. It works and I thank God! *D.H.*

Hip Pains – I've suffered with hip pains for years but never used any pain killer for good reasons. My husband bought a LLLT Pen for me and it removed all my pains by the third day. I use it almost every day to be safe now. *Vicky W*

Knee Problems – I've lived in pain with my left knee for decades. My right knee was replaced years ago and I was saving up \$50,000 to get a new left knee. A friend in the medical business said here Jim use this on your knee. I did and the pain seems to have left me for good, I think, I hope and pray. But of course I continue to use it daily. *Jim J*

Migraine Headaches – All my life, off and on, I have suffered with headaches and worst of all Migraine to the point of having to stay in a darkened room for days on end seeking relief. I have used both the Wand and the Disk behind my ears between the spinal column, and on the spinal column and normally within 30 minutes I have relief. It is rare that I have these headaches anymore but I keep my Pen and Disk handy, as they do work for me. *JS*

Eczemas – My 7 year old daughter has suffered all winter long with Eczemas. Her skin was red, itchy and many other problems pertaining to it. Nothing, absolutely nothing we have done helped her, that is, not until we got a treatment called "TLC-Disks". We used 2-Disks, strapped them on her two worst areas with Ace Bandages (AB) and put her to bed, in the morning she was 50% better, now each night we AB one on each foot, she's healed and no more itching. *N.A.*

Pains and sicknesses in our house – After years of most everyone in our family being sick and/or in pain for one thing and then another we used all of the TLC protocol, Dietary Supplements, Pens, Disks, the TLC Spectrabrite™ Clinical Unit, and now we are not sickly anymore, period! We recommend it to everyone. *Anthony and family*

Low Testosterone – At my age and that of many of my friends, we're not always able to perform our marital manly duties. A friend of mine told me use this TLC Disk. Put it in a small zip lock bag, put it in my underpants, and aim it face up at my privates all day if I wanted to, I did, and boy did it work. *Jim and others testify to this.*

DOCTORS: Countless testimonies have come back to us from individuals and even many doctors on how well our TLC Spectrabrite™ technologies have helped both the doctors, family members and their patients. To learn of more testimonies please go to our website: www.tlchealthsystems.com As we have time we'll continue to upload as many new testimonies as possible and we look forward to you and your patients sending us their testimonies as well.

OWNERS NOTES:

Treatment Comparison

Which One Would You Choose?

If there were a way you could own the TLC Spectrabrite™ Disk or Pen and use them as often as desired for pain relief for years to come (without side effects!) for a cash price of a month's supply of prescription painkillers or pain patches (with their stated potentially-*dangerous* side effects), which one would you choose?



Order TLC Health Systems Dietary Supplements presented above as needed to help buildup and sustain your total health. Our Essential Vitamins and Minerals shown are very unique, they are tasteless, children of all ages love them, and better yet these are Ionic, Fulvic and Humanic Minerals able to penetrate the bodies cell walls bringing life.